## Contents

Introduction
Storage and Stability
Kit Contents
Materials to Be Provided by User 3
Before Starting 3
Mag-Bind <sup>2</sup> Soil DNA Protocol
Mag-Bind <sup>?</sup> Soil DNA Clean-up Protocol
Troubleshooting Guide

#### Introduction

The Mag-Bind<sup>?</sup> Soil DNA Kit allows rapid and reliable isolation of high-quality total DNA from various soil samples. Up to 0.5 gram of soil samples can be processed in less than 60 minutes. The system combines the reversible nucleic acid-binding properties of MagSi<sup>?</sup> particles with the efficiency of HTR Reagents to eliminate humic acid, polysaccharides, phenolic compounds, and enzyme inhibitors from soil samples. Purified DNA is suitable for PCR, restriction digestion, and hybridization techniques. There are no organic extractions, thus reducing plastic waste and hands-on time to allow multiple samples to be processed in parallel.

## Storage and Stability

Most components of the Mag-Bind<sup>?</sup> Soil DNA Kit should be stored at 22°C-25°C. MagSi<sup>?</sup> Particles Solution should be stored at 4°C for long-term use. During shipment or storage in cool ambient conditions, precipitates may form in Buffer DS. Dissolve such deposits by warming the solution at 37°C and gently shaking.

## **Kit Contents**

Product Number	M5636-00	M5636-01	M5636-02
Purification	1 x 96	4 x 96 Preps	20 x 96 Preps
MagSi <sup>?</sup> Particles Solution	4.5 ml	18 ml	90 ml
Buffer SLX Mlus	80 ml	2 x 150 ml	6 x 250 ml
Buffer DS	10 ml	30 ml	150 ml
Buffer SP2	20 ml	80 ml	2 x 250 ml
cHTR Reagent	20 ml	80 ml	2 x 250 ml
Buffer XP5	120 ml	2 x 250 ml	10 x 250 ml
Buffer PHB	100 ml	2 x 200 ml	7 x 250 ml
SPM Wash Buffer	50 ml	4 x 50 ml	20 x 50 ml
Glass Beads	60 g	220 g	5 x 220 g
Elution Buffer	20 ml	80 ml	2 x 250 ml
Instruction Booklet	1	1	1

## **Before Starting**

Dilute SPM Wash Buffer with absolute ethanol as follows and store at room temperature.

M5636-00	Add 200 ml absolute (96%-100%) ethanol.
M5636-01	Add 200 ml (96%-100%) ethanol to each bottle.
M5636-02	Add 200 ml (96%-100%) ethanol to each bottle.

#### Materials to Be Provided by User

- Microcentrifuge capable of at least 14,000 x g
- Nuclease-free 1.5 ml or 2 ml microfuge tubes
- Water bath equilibrated to 65°C
- Absolute (96%-100%) ethanol

Please read the entire booklet to become familiar with Mag-Bind? Soil DNA Kit protocol.

# Mag-Bind<sup>?</sup> Soil DNA Kit Protocol

- Weigh 500 mg of glass beads in a 2 ml centrifuge tube, add 0.25-0.5 g soil sample. Add 0.7 ml Buffer SLX Mlus. Vortex at maximum speed for 5-10 minute to lyse samples. For the best result, A Mixer Mill, such as Fastprep-24<sup>?</sup>, Mixer Mill MM 300<sup>?</sup>, should be used.
- 2. Add 70 µl Buffer DS and vortex to mix.
- Incubate at 70°C for 10 min. Briefly vortex the tube once during the incubation. For some difficult lysis bacterial, Increase the temperature to 90°C.
- 4. Centrifuge at 13,000 x g for 5 min.
- 5. Transfer 500 μl the supernatant into a new 2 ml tube and add 170 μl Buffer SP2. Vortex to mix well.
- Add 170 µI cHTR Reagent to the sample. Mix throughly by vortexing for 10 seconds. Incubate on ice for 5 minutes. Centrifuge at 13,000 x g in a microcentrifuge for 5 minutes.
- Transfer 450 µl of the cleared supernatant to a new 1.5 ml tube or 2ml Deep Well Plate.
- Add 450 μl Buffer XP5 and 40 μl of MagSi<sup>?</sup> Particles to the sample. Shaking 60 seconds to mix well. Incubate at room temperature for 2 minutes.
- Place the tube or plate on a magnetic separation device suitable for 2 ml tube to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
- 10. Remove the tube containing the MagSi<sup>?</sup> Particles from the magnetic separation device.
- Add 500 μl Buffer XP5 into the tube or the plate. Resuspend MagSi<sup>?</sup> Particles pellet by vortexing.

3

```
4
```

- Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
  Add 800 µl Buffer PHB into the tube. Resuspend MagSi<sup>?</sup> Particles pellet by vortexing.
- Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
- 14. Add 800 μl of SPM Wash Buffer diluted with ethanol into the tube. Resuspend Mag-Bind<sup>?</sup> Particles pellet by vortexing.
- 15. Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
- Wash MagSi<sup>?</sup> Particles with SPM Wash Buffer one more time by repeating step 14-15.
- 17. After remove the supernatant, air dry the magnetic beads by invert the tube on a absorbent paper for 15 minutes. Remove any residue liquid from tube with pipettor.
- Add 50-100 ul Elution Buffer or water to the tube. Incubate the tube at Resuspend MagSi<sup>?</sup> Particles by vortexing. Incubate at 65°C for 10 minutes if maxium DNA yield is desired.
- Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles.
- 20. Transfer the cleared supernatant containing purified DNA to a new 1.5 ml tube.

## Mag-Bind<sup>?</sup> Soil DNA Clean-up Protocol

The following protocol is designed for clean-up the DNA purified with some inhouse purification methods or other commercial products.

- Transfer 200µl of crude DNA sample to a 1.5 ml tube. If the sample volume is less than 200ul, add Elution Buffer (supplied with this kit) or TE Buffer to bring the volume to 200µl.
- 2. Add 100μl of cHTR Reagent to the tube and Mix throughly by vortexing for 20 seconds.
- 3. Incubate at room temperature for 3 minutes.
- 4. Centrifuge at full speed (14,000-20,000 x g) for 2 minutes.
- 5. Transfer 250  $\mu$ l of the supernatant to a clean new 1.5 ml centrifuge tube, make sure not to disturb the pellet.
- Add 300 μl Buffer XP5 and 40 μl of MagSi<sup>?</sup> Particles to the sample. Invert to mix well. Incubate at room temperature for 2 minutes.
- Place the tube on a magnetic separation device suitable for 2 ml tube to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
- Remove the tube containing the MagSi<sup>?</sup> Particles from the magnetic separation device. Add 500 μl Buffer XP5 into the tube. Resuspend MagSi<sup>?</sup> Particles pellet by vortexing.
- Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant. Remove the tube containing the MagSi<sup>?</sup> Particles from the magnetic separation device.

5

6

- Add 800 μl Buffer PHB into the tube. Resuspend MagSi<sup>?</sup> Particles pellet by vortexing.
- Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
- 12. Add 800 μl of SPM Wash Buffer diluted with ethanol into the tube. Resuspend MagSi<sup>?</sup> Particles pellet by vortexing.
- Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles. Carefully remove and discard the cleared supernatant.
- 14. Wash MagSi<sup>?</sup> Particles with SPM Wash Buffer one more time by repeating step 12-13.
- 15. After remove the supernatant, air dry the MagSi<sup>?</sup> Particles by invert the tube on a absorbent paper for 15 minutes. Remove any residue liquid from tube with pipettor.
- Add 50-100 ul Elution Buffer or water to the tube. Incubate the tube at resuspend MagSi<sup>?</sup> Particles by vortexing. Incubate at 65°C for 10 minutes if maxium DNA yield is desired.
- 17. Place the tube onto a magnetic separation device to magnetize the MagSi<sup>?</sup> Particles.
- Transfer the cleared supernatant containing purified DNA to a new 1.5 ml tube.

## **Trouble Shooting Guide**

Problem	Cause	Suggestions	
Low DNA yield	Incomplete disruption of starting material.	Repeat the experiment with new sample, make sure the sample are completely interrupted and lysed.	
	Sample stored incorrectly.	Store the sample at -20°C.	
	Loss the MagSi <sup>?</sup> Particles during operation.	Carefully avoid remove the MagSi <sup>?</sup> Particles during aspiration.	
	DNA remains bound to MagSi <sup>?</sup> Particles.	Increase elution volume and incubate at 65°C for additional 10 min elution.	
	DNA washed off.	Dilute SPM Wash Buffer by adding appropriate volume of absolute ethanol prior to use (page 3).	
Problems in downstream applications		Wash Buffer must be at room temperature.	
	Ethanol carry-over.	Dry the MagSi <sup>?</sup> Particles before elution.	
	Add BSA to a final concentration of 0.1 μg/mL to the PCR mixture.	Add BSA to a final concentration of 0.1µg/mL to the PCR mixture.	
A260/280 ratio is low	Inefficient elimination of inhibitory compounds.	Repeat with a new sample, be sure to mix cHTR Reagent throughly before use.	
Eluted DNA still show color	Inefficient elimination of inhibitory compounds.	Start with a new sample, use extra cHTR or duplicate the cHTR extraction step.	

8

7